



AuPreP Random Primer DNA Labeling Mix System

for 25 labeling assays

Cat No. AUP-RPLM-25

The AuPreP Random Primer DNA Labeling Mix System is a specially designed premixed solution for the labeling of DNA with radiolabeled dCTP using random sequence oligonucleotides.

System Details:

The **AuPreP Random Primer Labeling Mix System** is based on the hybridization of oligonucleotides of all possible sequences to the denatured template DNA to be labeled. The complementary DNA strand is synthesized by a “klenow” fragment of DNA Polymerase I, using the random oligonucleotides as primers. By substituting a radiolabeled nucleotide for a non-radioactive equivalent in the reaction mixture, the newly synthesized complementary DNA is made radioactive.

The labeling mix system is a specially developed reaction mixture for enhanced convenience and performance. The reaction mixture contains random oligonucleotides, a Klenow fragment of DNA Polymerase I, dATP, dGTP, dTTP and a reaction buffer concentrate. The DNA labeling mix allows the labeling of the template DNA to a specific activity of 2×10^9 dpm/ μ g after only 10 minutes of incubation. This rapid labeling is accomplished with the use of the Klenow fragment, which lacks 5'-3' exonuclease activity, and by the use of nonamer primers giving more efficient priming from the template at 37°C. The labeling mix method enables the labeling of small amounts of DNA (10-20ng), such as restriction fragments isolated from gels. Fragments can be labeled directly in low melting temperature agarose gel slices. The labeled probes are used in various hybridization techniques, such as Southern and Northern blots, in-situ hybridization and screening of gene libraries.

System Components:

One vial containing 100 μ l DNA labeling mixture for 25 labeling assays.

Storage and Stability: The premixed Solution should be stored at -20°C. Avoid repeated changes in the solution temperature.

System Protocol

Add 10-25ng of template DNA to be labeled to a specific activity of 2×10^9 dpm/ μ g. The mixture is designed for use with (alpha³² P) dCTP with a specific activity of 3000 ci/mmol.

Standard Labeling Procedure

Add 10-25ng of template DNA to be labeled and sterile double distilled water to a final volume of 11 µg in a microcentrifuge tube.

Denature the DNA sample by heating to 95-100°C for 5 minutes, and then chill quickly on ice for 5 minutes.

Mix on ice:

- 11 µl (10-25ng) of denatured DNA.
- 4 µl of Labeling Mix Solution.
- 5 µl of (alpha³²P) dCTP (3000 ci/mmol).

Incubate at 37°C for 10 minutes.

Stop the reaction by adding 2µl of 0.2M EDTA (pH=8), or by heating to 65°C for 7 minutes.

For use in hybridization, denature the labeled DNA by heating to 95°C for 5 minutes and then cool on ice.

Important Notes:

- A probe with specific activity above 1×10^9 dpm/µg can be obtained after 3 minutes of incubation. The maximum probe specific activity is obtained after 10-20 minutes of incubation.
- Less than 10-20ng DNA can also be labeled, but the maximal incorporation may be achieved only after 30-60 minutes.
- DNA fragment in low melting temperature agarose can be used directly in the reaction without the removal of the agarose (see Appendix I).
- Removal of unincorporated nucleotides is not necessary, but sometimes it is desirable for reducing the background during hybridization. When required, the probes can be purified by gel filtration (Sephadex G-50) or by ethanol precipitation.
- Increasing the amount of template DNA above 25ng will reduce the specific activity of the labeled DNA.

Appendix I: Labeling of DNA fragments in low melting temperature agarose

- I After agarose gel electrophoresis, cut out a slice of the gel containing the target DNA fragment.
- II Add 3ml distilled sterilized water for each gram of the gel slice.
- III Place in boiling water for 10 minutes to melt the gel and to denature the DNA. Keep at 37°C until required.
- IV The DNA solution can be used directly in the reaction as template DNA (proceed from Step 1 of the standard labeling procedure).

<u>Other AuPreP™ DNA/RNA Kits</u>	<u>Other Related Products</u>
AuPreP™ Plasmid Maxi Kit	AuPreP Oligos (High Affinity Purified Oligo synthesis available in different scales, purifications & modifications)
AuPreP™ Plasmid Midi Kit	AuPreP TaQ DNA Polymerase (Ultrapure, Ultra-stable & Ultra-sensitive Taq DNA Polymerase)
AuPreP™ SPIN™ SPIN Miniprep Kit	AuPreP Hotstart TaQ DNA Polymerase (Robust Polymerase for Hotstart PCR assays)
AuPreP™ Blood Genomic DNA Maxi	AuPreP Super Fidelity TaQ DNA Polymerase (High fidelity Polymerase produces blunt ended amplicons upto 5Kb)
AuPreP™ Blood Genomic DNA Extraction Midi Kit	PCR Doctor - (PCR enhancer for AuPreP Hotstart Taq or Super Fidelity Taq especially designed for GC/AT/Dirty/Difficult Templates)
AuPreP™ GEN^{ht} DNA Extraction Kit	AuPreP Longjump Polymerase (Robust Long Polymerase for templates > 4kb to 18kb+ for challenging PCRs)
AuPreP™ DNA easy Plant Maxi kit	AuPreP Red PCR Master Mix (2x Master mix with Red Dye without Enhancer)
AuPreP™ DNA easy Plant Mini Kit	AuPreP DIAMOND MASTER-MIX (2x Mastermix with PCR Enhancer & Stabilizer without tracking dyes)
AuPreP™ PCR Purification Kit	AuPreP DIAMOND DOUBLE DYE MASTERMIX (2x Mastermix with PCR Enhancer, Stabilizer & tracking dyes)
AuPreP™ Plant RNA Maxi Kit	AuPreP DNA Extraction System (A fast Reagent for pure genomic DNA isolation for down stream applications)
AuPreP™ Plasmid Maxi Kit	AuPreP RNA Extraction System (for Purest & High Quality RNA extraction with simple cost effective protocol)
AuPreP™ RNA Easy Midi Kit	AuPreP Gold cDNA Synthesis Kit (Highly Cost effective cDNA Synthesis Kit using RT with reduce Rnase H activity)
AuPreP™ RNA™ Mini Kit	AuPreP Gold RT-PCR Combo Kit (2 step RT-PCR protocol with tracking Dye)
AuPreP™ RNV™ Viral RNA Extraction Miniprep Kit	AuPreP Extra Mile First Strand cDNA System (Premium cDNA Synthesis Kit using RT with point mutant Rnase H minus activity)
	Novascript III RNase H⁻ RT (Premium Ultra-stable Rnase H minus RT for long high quality cDNA construction)
	Novascript III single step RT-PCR System (Premium 1step RT-PCR system using Novascript & AuPreP Hotstart DNA Polymerase)
	AuPreP Random Primer labeling Mix System (Premixed solution for the labeling of DNA with radiolabeled dCTP using random sequence oligonucleotides)

References

1. Feinberg, A.P. and Vogelstein, B., *Anal. Biochem.*, **132**:6-13, 1983
2. Feinberg, A.P. and Vogelstein, B., *Anal. Biochem.*, **137**:266-267, 1984
3. Grunstein, M. And Hogness, D.S., *Proc. Nati. Acad. Sci. USA*, **72**:3961-3965, 1975
4. Clark, J.M., Joyce, C.M. and Beardsley, G.P., *J. Mol. Biol.*, **198**:123-127, 1987