



## **AuPreP™ DNA easy Plant Maxi kit**

**Cat. # PLX68-163LTD**

Isolation of genomic DNA from 1 g plant material.

### **Kit Contents:**

PX1 Buffer	90 ml
PX2 Buffer	30 ml
PX3 Buffer	110 ml
RNase A	85 mg
WS Buffer	50 ml
Plant Genomic DNA Maxi Column	20
Shearing tube	20
15ml Collection tube	40
Protocol	1

### **Notes:**

- All centrifugation should be done at room temperature with a swing-bucket centrifuge.
- Preheat a water bath to 65°C.
- Preheat TE or ddH<sub>2</sub>O to 65°C for DNA elution.
- PX1 Buffer and PX3 Buffer may form a precipitate, warm at 65°C to redissolve.
- Add 850 µl of ddH<sub>2</sub>O to the RNase A powder tube, vortex to dissolve and store at 4°C.

### **Protocol:**

- 1. Grind 1 g (or less) plant sample under liquid nitrogen to a fine powder and transfer to a new tube.**  
Do not allow the sample to thaw, and continue immediately to step 2.
- 2. Add 4 ml of PX1 Buffer and 40 µl of RNase A solution (100 mg/ml) to the tissue powder and vortex vigorously, then incubate the mixture at 65°C for 10 minutes.**  
Do not mix PX1 Buffer and RNase A prior to use. Invert 2-3 times during 65°C incubation.
- 3. Add 1.3 ml of PX2 Buffer to the lysate, vortex, and incubate on ice for 5 minutes.**
- 4. Apply lysate to the Shearing tube sitting in a Collection tube and centrifuge at full speed (about 3000 rpm or 2500 x g) for 2 minutes. Transfer flow-through sample from the Collection tube to a new tube (not provided).**  
Avoid pipetting any debris or pellet in the collection tube.

To order contact:

306, Aggarwal City Mall, Opposite M2K Pitampura,  
Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222  
Mobile: +91-9810521400 Fax: +91-11-42208444  
Email: [customerservice@atzlabs.com](mailto:customerservice@atzlabs.com)  
Web: [www.atzlabs.com](http://www.atzlabs.com)



- 5. Add 0.5 volume of PX3 Buffer and 1 volume of 96-100% ethanol to the clear lysate and mix by pipetting.**  
For example: If 4.5 ml clear lysate collected, add 2.25 ml PX3 Buffer and 4.5 ml ethanol.
- 6. Apply 6.5 ml of the ethanol added sample (including any precipitate) from step 5 to a Plant Genomic DNA Maxi Column sitting in a Collection tube, close the cap, centrifuge at full speed for 3 minutes, and discard the filtrate.**  
If the solution remains above the membrane, centrifuge again.
- 7. Repeat step 6 for rest of the sample.**
- 8. Wash the column twice with 5 ml of WS Buffer by centrifuging at full speed for 3 minutes and discard the filtrate.**  
Add 200 ml of ethanol (96-100%) to the WS Buffer bottle when first open the bottle.
- 9. Centrifuge at full speed for 5 minutes to remove traces of WS Buffer.**
- 10. Transfer the column to a new 15 ml tube (not provided), add 2 ml of 65°C TE or ddH<sub>2</sub>O, and centrifuge at full speed for 5 minutes to elute DNA.**
- 11. Store DNA at -20°C.**

**To order contact:**

306, Aggarwal City Mall, Opposite M2K Pitampura,  
Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222  
Mobile: +91-9810521400 Fax: +91-11-42208444  
Email: [customerservice@atzlabs.com](mailto:customerservice@atzlabs.com)  
Web: [www.atzlabs.com](http://www.atzlabs.com)