



8. Precautions

- 1) This product is suitable for fresh milk and pure milk.
- 2) The milk sample should store in cold, the storage time should be no more than 3 days.
- 3) Milk samples should be fully liquid without any agglomeration and deposition. Strongly recommend to mix sample evenly to ensure the result is accurately.
- 4) The test strips and microwells can be used only once, do not use test strips out of expiry date.
- 5) Do not touch the white membrane surface in the middle of test strip, avoid sunlight and fan blowing directly.
- 6) For sample with high fat content, Chromatography speed is slow, recommend to extend the second step incubation for another 1min.
- 7) Use the test dipsticks testing again for positive results.
- 8) Please contact the supplier for any questions.

9. Storage and expiry date

Storage: Store at 2-8 °C in dark, sealed, dry place, no frozen.

Expiry date: 12 months; date of production is on box.

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Road No. 44, Pitampura, Delhi – 110034, India

Mobile: +91-98105-21400, Tel: +91-11-42208000, 8111, 8222, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com, www.atzlabs.com ; www.lifetechindia.com