

PathoScreen® Kit for PVY-n

Test Label: Alkaline Phosphatase

Test Format: Triple Antibody Sandwich (TAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Monoclonal

Intended Use:

This test is used for detection of isolates of n strain group of PVY. This test also detects the ntn strains that cause the tuber necrotic ringspot disease (PTNRD) of potatoes.

Includes:

- Antibody-coated 96-well microtiter plates
- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- Nonfat dried milk
- PNP substrate tablets, 5 mg each
- PNP substrate buffer, 5X concentrate
- Positive control
- PBST buffer, 20X concentrate, 50 ml
- Tween-20
- [User Guide](#)
- [Tech Note](#)



PathoScreen® Kit, Potato virus Y -necrotic (PVY -n strain),PSA 26001/0096

PathoScreen® Kit, Potato virus Y -necrotic (PVY -n strain),PSA 26001/0288

PathoScreen Kit, Potato virus Y -necrotic (PVY -n strain),PSA 26001/0480

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com