

MPR-20	<u>ELISA Microplate Reader (automated, programmable) with 4 filters (no printer or software)</u>
MPR-20-6U	<u>ELISA Microplate Reader (automated, programmable) with 6 filters (UV range)</u>
MPR-20-6V	<u>ELISA Microplate Reader (automated, programmable) with 6 filters (Visible range)</u>
MPR-30	<u>ELISA Microplate Reader (automated, programmable) with 4 filters & ELISA reader software</u>
MPR-30-6U	<u>ELISA Microplate Reader (automated, programmable) with 6 filters (UV range) & ELISA reader software</u>
MPR-30-6V	<u>ELISA Microplate Reader (automated, programmable) with 6 filters & ELISA reader software</u>
MPR-30-6VP	<u>ELISA Microplate Reader (automated, programmable) with 6 filters (Visible range), on-board printer & ELISA reader software</u>
MPR-30-P	<u>ELISA Microplate Reader (automated, programmable) with 4 filters and on-board printer & ELISA reader software</u>
MPR-404	<u>ELISA Microplate Reader (compact, 8 channel, automated, programmable) with 4 filters & ELISA reader software</u>
MPR-406	<u>ELISA Microplate Reader (compact, 8 channel, automated, programmable) with 6 filters & ELISA reader software</u>
MPR-504	<u>ELISA Microplate Reader (compact, 8 channel, automated, programmable) with 4 filters & ELISA reader software with on-board printer</u>
MPR-506	<u>ELISA Microplate Reader (compact, 8 channel, automated, programmable) with 6 filters & ELISA reader software with on-board printer</u>
MPW-30	<u>ELISA Microplate Washer Automated, programmable</u>
MPW-60	<u>ELISA Microplate Washer, Manual</u>
MSR-10	<u>ELISA Microwell Strip Reader with built-in dot matrix printer</u>