

## Single Channel Mini-Pipettes (Fixed-Volume)

### Single Channel Mini-Pipettes

#### Single Channel Mini-pipette



**Sizes: 5, 10, 20, 25, 50, 100 ul**

- Cutest pipettes available
- Minimum size of 130 mm-no stress
- Excellent accuracy and precision
- Two stage plunger for last drop dispensing
- Smart Tip cone design allows perfect fit for 20 ul or 200 ul
- Autoclavable and Color coded
- Available in 5, 10, 20, 25, 50, 100 ul sizes
- Most economical and affordable

Catalog#	Product Description
SCMP-5	Single channel mini-pipette (130 mm) 5 ul, autoclavable
SCMP-10	Single channel mini-pipette (130 mm) 10 ul, autoclavable
SCMP-20	Single channel mini-pipette (130 mm) 20 ul, autoclavable
SCMP-25	Single channel mini-pipette (130 mm) 25 ul, autoclavable
SCMP-50	Single channel mini-pipette (130 mm) 50 ul, autoclavable
SCMP-100	Single channel mini-pipette (130 mm) 100 ul, autoclavable
SCMP-Set5	Single channel mini-pipette (130 mm) set of 5 pipettes(5, 10, 20, 50, 100 ul), autoclavable (note: 1 free pipette of 25 ul included)