

PathoScreen[®] Kit for PNRSV

Test Label: Alkaline Phosphatase

Test Format: Triple Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Monoclonal

Intended Use:

This ELISA is a qualitative serological assay for detection of *Prunus necrotic ringspot virus* (PNRSV) in stone fruits.

Symptoms of PNRSV can include stunting, chlorosis, and leaf deformities. PNRSV is also the virus most associated with *Rose mosaic disease*.

Includes:

- Antibody-coated 96-well microtiter plates
- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- ECI buffer, 5X concentrate
- PNP substrate tablets, 5 mg each
- PNP substrate buffer, 5X concentrate
- Positive control
- PBST buffer, 20X concentrate, 50 ml
- Tween-20
- General extract buffer, GEB
- [User Guide](#)



PathoScreen[®] Kit, Prunus necrotic ringspot virus (PNRSV), PSA 30500/0096

PathoScreen[®] Kit, Prunus necrotic ringspot virus (PNRSV), PSA 30500/0288

PathoScreen[®] Kit, Prunus necrotic ringspot virus (PNRSV), PSA 30500/0480

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com