



BIOLOGICAL INDUSTRIES

PET media bottles (pre-sterilized)

These square PET media bottles offer the clarity of glass with the durability of plastic and a space saving design. Improved features include a wide mouth opening, thicker wall dimensions and a leak proof closure. Bottles and caps are pre-sterilized and tested for pH stability, temperature durability, cloning efficiency and cytotoxicity.

Product Name	Size
PET Plastic (Media) Bottles (pre-sterilized)	20 ml (612/carton)
	100 ml (96/carton)
	500ml (48/carton)

To order contact:
306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222
Mobile: +91-9810521400 Fax: +91-11-42208444
Email: customerservice@atzlabs.com
Web: www.atzlabs.com