

To order contact:
306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222
Mobile: +91-9810521400 Fax: +91-11-42208444
Email: customerservice@atzlabs.com
Web: www.atzlabs.com



Contains: (25) Tubes RNAsable, (1) Tube Sterile Water and Quick Reference Protocol. Enables long-term stabilization of RNA at room temperature. Easy sample recovery by simply adding water.



Contains: (2) RNAsable 96-Well Plates, (2) Tubes Sterile Water and Quick Reference Protocol. Enables long-term stabilization of RNA at room temperature. Easy sample recovery by simply adding water.