

## ELISA Reagent Set for Apple mosaic virus (ApMV)

**Test Label:** Alkaline Phosphatase

**Test Format:** Double Antibody Sandwich (DAS) ELISA

**Capture Reagent:** Polyclonal

**Detection Reagent:** Polyclonal

### Intended Use:

Agdia's ApMV DAS-ELISA test system is a qualitative serological assay for detection of *Apple mosaic virus* (ApMV) in plant tissue with or without external disease symptoms.

The test can be used by growers, plant breeders, disease clinics, researchers, commercial lab services, regulatory agencies, and anyone testing for ApMV infection.

### Includes:

- Capture antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid\* (\*strip in single plate set)
- [User Guide](#)



Reagent Set for Apple mosaic virus (ApMV), 96 testwells, SRA 30700/0096

Reagent Set for Apple mosaic virus (ApMV), 500 testwells, SRA 30700/0500

Reagent Set for Apple mosaic virus (ApMV), 1000 testwells, SRA 30700/1000

Reagent Set for Apple mosaic virus (ApMV), 5000 testwells, SRA 30700/5000

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi - 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)