

## ELISA Reagent Set for Carnation etched ring virus (CERV)

**Test Label:** Alkaline Phosphatase

**Test Format:** Triple Antibody Sandwich (TAS) ELISA

**Capture Reagent:** Polyclonal

**Detection Reagent:** Monoclonal with Anti-animal conjugate

### Intended Use:

This ELISA test is a qualitative serological assay used for detection of *Carnation etched ring virus* (CERV) in Carnation and Dianthus.

### Includes:

- Capture antibody
- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- 96-well microtiter plates, solid\* (\*strip in single plate set)
- [User Guide](#)



Reagent Set for Carnation etched ring virus (CERV), 96 testwells, SRA 34100/0096

Reagent Set for Carnation etched ring virus (CERV), 500 testwells, SRA 34100/0500

Reagent Set for Carnation etched ring virus (CERV), 1000 testwells, SRA 34100/1000

Reagent Set for Carnation etched ring virus (CERV), 5000 testwells, SRA 34100/5000

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)